


cardiolife connect

BECAUSE IT MATTERS
SMARTPHONE APP
GPS TRACKING
24/7 AED SERVICE



CARDIAC SAFETY

Nihon Kohden's cardiolife connect solution provides full-service AED management, including complete monitoring of the cardiolife AED-3100 and other connected AEDs. Not only do you gain insight into the AED's current status, but also its exact location (GPS) and whether it has been used to rescue a patient.

The entire platform is customizable to individual requirements in terms of notifications and integration.

24/7 SURVEILLANCE AND READINESS

cardiolife connect is used for 24/7 surveillance of the AED. Nihon Kohden provides monitoring of the AED to verify that the device is operational, such as the shelf life of the disposable pads or the battery level.

You can customize whether you want to receive alerts for daily and monthly self-tests, or when specific actions occur, such as when the AED is turned on or when the AED leaves a certain location.

cardiolife connect allows you to select different notifications.

You can also set whether you want to receive email or SMS notifications about these self-selected events.




E-mail



SMS





Messenger



DASHBOARD AND APPS

cardiolife connect's technology consists of a communication module and a web-accessible dashboard that provides a bird's-eye view of all AEDs, allowing the operator the ability to monitor the service level and readiness of all deployed units at all times. This single, unified view can also be accessed via the CardiLink® apps for iOS and Android smartphones.




The app intuitively facilitates the data visualisation. Registration is made the easiest way possible via the CardiLink® 2.0 app, available for both Android and iOS devices.





GPS TRACKING AND GEO FENCING

cardiolife connect enables real GPS tracking. Exact information about the AED's location is displayed on a digital map.




Motion activity details are available in the movement list. The GEO fencing function triggers an alarm if the AED leaves a certain defined radius.



INTERFACE AND FLEET MANAGEMENT

Due to the interface cardioline connect offers a unique integration which is able to transfer data to mobile first responder applications.


Our fleet management capabilities connect the sometimes disparate technologies deployed by health and safety device manufacturers on one networked platform – mitigating complications arising from local device deployments.



PRODUCT COMPATIBILITY

This solution offers full hardware compatibility. Nihon Kohden's current cardioline AED-3100 defibrillator can be equipped with the cardiolineconnect solution.

cardioline **AED**
AED-3100





Improving Healthcare with Advanced Technology

Since its foundation back in 1951, Nihon Kohden's mission has been to improve the quality of life with advanced technology. We provide solutions for diagnosis, critical care, clinical information, and in vitro diagnostics – and we are dedicated to collaborate with you to meet the challenges of healthcare today and tomorrow.

Visit www.nihonkohden.com to find out more.

CardiLink GmbH is the technology partner for
Nihon Kohden's cardiolifeconnect solution.

NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, 61191 Rosbach, Germany
Phone: +49 6003 827 0, Fax: +49 6003 827 599
Internet: <https://eu.nihonkohden.com>, E-mail: info@nke.de

CardiLink GmbH
Gebhardtstraße 28, 6. OG OST, 90762 Fürth, Germany
Phone: +49 911 733 050 44
Internet: www.cardi-link.com



This document may be revised or replaced by
Nihon Kohden at any time without notice.