

AC Adapter

SC-101V

General

The SC-101V AC adapter is for the EMS-1052 defibrillator and SB-101V battery charger. An optional power cord is used together with the adapter.

Safety Information

⚠ WARNING A warning alerts the user to possible injury or death associated with the use or misuse of the instrument.

⚠ CAUTION A caution alerts the user to possible injury or problems with the instrument associated with its use or misuse such as instrument malfunction, instrument failure, damage to the instrument, or damage to other property.

Pay attention to all safety information in the Operator's Manual or Installation Guide.

⚠ WARNING

Do not disassemble or modify the AC adapter. It might cause skin burn, fire, electrical shock or injury.

⚠ WARNING

- When connecting the EMS-1052 defibrillator to an AC outlet in an ambulance using the AC adapter, make sure that all AC outlets are securely grounded and have no defects. Otherwise noise may interfere with the power or the operator may receive electrical shock.
- The wiring and grounding differ depending on the ambulance. Therefore, Nihon Kohden cannot provide instructions for connecting the defibrillator. You must have a qualified electrician install the AC outlet.

⚠ CAUTION

Connect the AC adapter only to specified instruments and follow the specified procedure. Failure to follow this instruction may result in electrical shock or injury to the operator due to leakage current, or may cause fire or damage to the AC adapter.

This Safety and Performance Information is an extract from the general and safety information sections of the most recent edition of Operator's Manual or Installation Guide. Therefore, the contents of your Operator's Manual or Installation Guide may differ from those of this Safety and Performance Information. For detailed operating procedures, follow the instructions of your Operator's Manual or Installation Guide.

 Manufacturer
NIHON KOHDEN CORPORATION
1-31-4 Nishiochiai, Shinjuku-ku,
Tokyo 161-8560, Japan
Phone +81 3-5996-8041
<https://www.nihonkohden.com/>

EC REP European Representative
NIHON KOHDEN EUROPE GmbH
Raiffeisenstrasse 10, D-61191 Rosbach, Germany
Phone +49 6003-827-0 Fax +49 6003-827-599