

MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS23X
	2023-11-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



		MEK -6400 -6410 -6420	MEK -6500 -6510	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305	MEK -6450	MEK -6550
WBC	10 ³ /μL	9.3	9.3	9.6	9.4	9.3	9.69	9.27	9.3	9.3
RBC	10 ⁶ /μL	4.70	4.70	4.70	4.70	4.69	4.72	4.69	4.70	4.70
HGB	g/dL	13.9	14.0	13.6	13.7	13.5	13.77	13.76	13.9	14.0
HCT	%	39.6	39.6	41.9	40.7	39.5	46.4	46.1	39.6	39.6
MCV	fL	84.2	84.2	89.1	86.6	84.3	98.2	98.4	84.2	84.2
MCH	pg	29.6	29.8	28.9	29.1	28.8	29.2	29.3	29.6	29.8
MCHC	g/dL	35.1	35.4	32.5	33.7	34.2	29.7	29.8	35.1	35.4
PLT	10 ³ /μL	253	253	234	237	236	259.8	269.5	253	253
MPV	fL	7.1	9.0	7.1	7.8	7.1	9.5	9.5	7.1	9.0
RDW-CV *	%	15.8	15.8	14.7	16.0	15.0	15.7	15.6	15.8	15.8

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	170
FL	85
SD	113
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901780H