


MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS23Y
	2023-12-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



<https://www.nihonkohden.com/spi/index.html>



 **NIHON KOHDEN CORPORATION**
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

		MEK -6400 -6410 -6420	MEK -6500 -6510	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305	MEK -6450	MEK -6550
WBC	10 ³ /μL	9.3	9.3	9.4	9.3	9.1	9.50	9.27	9.3	9.3
RBC	10 ⁶ /μL	4.58	4.58	4.55	4.52	4.55	4.60	4.57	4.58	4.58
HGB	g/dL	13.6	13.7	13.4	13.4	13.3	13.59	13.57	13.6	13.7
HCT	%	39.4	39.4	41.4	40.2	39.2	46.0	45.7	39.4	39.4
MCV	fL	86.1	86.1	91.0	88.9	86.1	100.1	99.9	86.1	86.1
MCH	pg	29.7	29.9	29.5	29.6	29.2	29.5	29.7	29.7	29.9
MCHC	g/dL	34.5	34.8	32.4	33.3	33.9	29.5	29.7	34.5	34.8
PLT	10 ³ /μL	247	247	232	234	234	257.7	262.2	247	247
MPV	fL	6.7	8.8	6.7	7.6	6.7	9.4	9.3	6.7	8.8
RDW-CV *	%	16.1	16.1	15.0	16.3	15.3	16.1	16.1	16.1	16.1

MEK	
	-7222
	-7300
	-8222
	-9100
	-9200
FS	168
FL	84
SD	106
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



0694-901780I