


MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS23Z
	2024-01-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



		MEK -6400 -6410 -6420	MEK -6500 -6510	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305	MEK -6450	MEK -6550
WBC	10 ³ /μL	9.0	9.0	9.1	9.0	8.9	9.19	9.02	9.0	9.0
RBC	10 ⁶ /μL	4.68	4.68	4.65	4.63	4.65	4.70	4.67	4.68	4.68
HGB	g/dL	13.8	14.0	13.7	13.7	13.5	13.83	13.79	13.8	14.0
HCT	%	40.4	40.4	41.9	40.7	39.7	46.3	46.0	40.4	40.4
MCV	fL	86.4	86.4	90.2	88.0	85.3	98.5	98.6	86.4	86.4
MCH	pg	29.5	29.9	29.5	29.6	29.0	29.4	29.5	29.5	29.9
MCHC	g/dL	34.2	34.7	32.7	33.7	34.0	29.9	30.0	34.2	34.7
PLT	10 ³ /μL	253	253	239	242	241	267.9	271.8	253	253
MPV	fL	7.2	9.1	7.2	7.8	7.2	9.6	9.5	7.2	9.1
RDW-CV *	%	16.6	16.6	15.2	16.7	15.6	16.5	16.4	16.6	16.6

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	167
FL	85
SD	107
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901780J