


**MK-3CH** HEMATOLOGY CONTROL+CRP 2.5mL

**ASSAY SHEET**

<b>LOT</b>	<b>MK3C238H</b>
	<b>2023-11-05</b>
<b>OPEN</b>	<b>14DAYS</b>

**MEK-1303 Barcode Entry**

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 <sup>3</sup> /μL	19.48	± 2.20
RBC	10 <sup>6</sup> /μL	5.64	± 0.25
HGB	g/dL	18.73	± 0.70
HCT	%	60.2	± 4.0
MCV	fL	106.7	± 7.0
MCH	pg	33.2	± 2.8
MCHC	g/dL	31.1	± 3.0
PLT	10 <sup>3</sup> /μL	567.3	± 70.0
MPV	fL	9.1	± 3.0
RDW-CV	%	16.0	± 3.0
LY%	%	16.90	± 10.00
MO%	%	0.82	± 5.00
GR%	%	82.28	± 10.00
LY	10 <sup>3</sup> /μL	3.29	± 1.80
MO	10 <sup>3</sup> /μL	0.16	± 0.70
GR	10 <sup>3</sup> /μL	16.03	± 2.00
CRP	Precision	≤ 10.0%	
		(41.10-76.30 mg/L)	

Manual Registration	
Lot	238
Level	High
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**

 1-31-4 NISHIOCHIAI, SHINJUKU-KU  
 TOKYO 161-8560 JAPAN

**0694-901767C**