


MK-3CH HEMATOLOGY CONTROL+CRP 2.5mL

ASSAY SHEET

LOT	MK3C23YH
	2024-02-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 ³ /μL	19.82	± 2.20
RBC	10 ⁶ /μL	5.74	± 0.25
HGB	g/dL	18.81	± 0.70
HCT	%	60.6	± 4.0
MCV	fL	105.5	± 7.0
MCH	pg	32.8	± 2.8
MCHC	g/dL	31.0	± 3.0
PLT	10 ³ /μL	582.9	± 70.0
MPV	fL	8.9	± 3.0
RDW-CV	%	16.5	± 3.0
LY%	%	14.75	± 10.00
MO%	%	1.88	± 5.00
GR%	%	83.37	± 10.00
LY	10 ³ /μL	2.92	± 1.80
MO	10 ³ /μL	0.37	± 0.70
GR	10 ³ /μL	16.53	± 2.00
CRP	Precision	≤ 10.0%	
		(43.40-80.60 mg/L)	

Manual Registration	
Lot	23Y
Level	High
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN

0694-901767D