


**MK-3CL** HEMATOLOGY CONTROL+CRP 2.5mL

**ASSAY SHEET**

<b>LOT</b>	<b>MK3C238L</b>
	<b>2023-11-05</b>
<b>OPEN</b>	<b>14DAYS</b>

**MEK-1303 Barcode Entry**

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 <sup>3</sup> /μL	2.01	± 0.60
RBC	10 <sup>6</sup> /μL	2.19	± 0.15
HGB	g/dL	5.72	± 0.40
HCT	%	18.9	± 2.0
MCV	fL	86.3	± 9.0
MCH	pg	26.1	± 2.4
MCHC	g/dL	30.3	± 3.0
PLT	10 <sup>3</sup> /μL	70.2	± 30.0
MPV	fL	8.6	± 3.0
RDW-CV	%	15.8	± 3.0
LY%	%	62.43	± 10.00
MO%	%	5.07	± 5.00
GR%	%	32.50	± 10.00
LY	10 <sup>3</sup> /μL	1.25	± 0.50
MO	10 <sup>3</sup> /μL	0.10	± 0.20
GR	10 <sup>3</sup> /μL	0.66	± 0.50
CRP	Precision	≤ 10.0%	
		(1.10-6.30 mg/L)	

Manual Registration	
Lot	238
Level	Low
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**

 1-31-4 NISHIOCHIAI, SHINJUKU-KU  
 TOKYO 161-8560 JAPAN

**0694-901765C**