


**MK-3CN** HEMATOLOGY CONTROL+CRP 2.5mL


**ASSAY SHEET**

<b>LOT</b>	<b>MK3C235N</b>
	<b>2023-08-05</b>
<b>OPEN</b>	<b>14DAYS</b>

**MEK-1303 Barcode Entry**

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 <sup>3</sup> /μL	7.62	± 0.80
RBC	10 <sup>6</sup> /μL	4.62	± 0.20
HGB	g/dL	13.79	± 0.50
HCT	%	45.0	± 3.0
MCV	fL	97.3	± 7.0
MCH	pg	29.8	± 2.4
MCHC	g/dL	30.6	± 3.0
PLT	10 <sup>3</sup> /μL	306.2	± 50.0
MPV	fL	9.0	± 2.0
RDW-CV	%	16.4	± 3.0
LY%	%	33.05	± 10.00
MO%	%	3.32	± 5.00
GR%	%	63.63	± 10.00
LY	10 <sup>3</sup> /μL	2.52	± 1.00
MO	10 <sup>3</sup> /μL	0.25	± 0.40
GR	10 <sup>3</sup> /μL	4.85	± 1.00
CRP	Precision	≤ 10.0% (17.30-32.10 mg/L)	

Manual Registration	
Lot	<b>235</b>
Level	<b>Normal</b>
Control	<b>MK-3C</b>


<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**

 1-31-4 NISHIOCHIAI, SHINJUKU-KU  
 TOKYO 161-8560 JAPAN


0694-901766B