


MK-RE1 HEMATOLOGY CONTROL
for reticulocyte **3.0mL**

ASSAY SHEET

LOT	MKRE239-1
	2023-11-05
OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	1.22	±1.00
RET	10⁶/μL	0.0433	±0.0355

2D Barcode Entry



• If RBC value is under 3.00 [$\times 10^6/\mu\text{L}$], please refrain from using.
It may have deteriorated or hemolyzed, which affects the measurement accuracy.



<https://www.nihonkohden.com/spi/index.html>

Manual Registration						
Lot	239					
Level	1					
Expiration	11	M	5	D	2023	Y
Control	MKRE					



NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN