


MK-RE2
HEMATOLOGY CONTROL
 for reticulocyte

3.0mL

ASSAY SHEET

| | |
|---|------------|
| LOT | MKRE239-2 |
|  | 2023-11-05 |
| OPEN | 14 DAYS |

| | | | |
|-------------|--------------------------|-----------------|----------------|
| | | MEK-9200 | RANGE |
| RET% | % | 3.42 | ±3.25 |
| RET | 10⁶/μL | 0.1204 | ±0.1145 |

2D Barcode Entry


• If RBC value is under 3.00 [$\times 10^6/\mu\text{L}$], please refrain from using.
 It may have deteriorated or hemolyzed, which affects the measurement accuracy.


<https://www.nihonkohden.com/spi/index.html>
Manual Registration

| | | | | | | |
|------------|------|---|---|---|------|---|
| Lot | 239 | | | | | |
| Level | 2 | | | | | |
| Expiration | 11 | M | 5 | D | 2023 | Y |
| Control | MKRE | | | | | |


 0123



NIHON KOHDEN CORPORATION
 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN