

**MEK-CAL** HEMATOLOGY CALIBRATOR 2.0mL

**ASSAY SHEET**

LOT	PLUS241
	2024-02-05
OPEN	7DAYS

MEK-9100,9200  
2D Barcode Entry



MEK-1301,1302,1303,1305  
Barcode Entry



		MEK -6400 -6410 -6420	MEK -6500 -6510	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305	MEK -6450	MEK -6550
WBC	10 <sup>3</sup> /μL	9.1	9.1	9.3	9.2	9.1	9.41	9.19	9.1	9.1
RBC	10 <sup>6</sup> /μL	4.63	4.63	4.61	4.60	4.61	4.65	4.64	4.63	4.63
HGB	g/dL	13.6	13.7	13.4	13.4	13.3	13.61	13.60	13.6	13.7
HCT	%	39.6	39.6	41.2	40.1	39.0	46.0	45.9	39.6	39.6
MCV	fL	85.6	85.6	89.3	87.1	84.5	98.9	98.9	85.6	85.6
MCH	pg	29.4	29.6	29.1	29.1	28.9	29.3	29.3	29.4	29.6
MCHC	g/dL	34.3	34.6	32.5	33.4	34.1	29.6	29.6	34.3	34.6
PLT	10 <sup>3</sup> /μL	255	255	249	253	252	263.3	270.6	255	255
MPV	fL	7.0	8.9	7.0	7.7	7.0	9.3	9.3	7.0	8.9
RDW-CV *	%	15.9	15.9	14.8	16.1	15.1	15.9	15.8	15.9	15.9

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	166
FL	89
SD	108
FS THR.	100

\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN



0694-901780K