


MEK-5DL HEMATOLOGY CONTROL
for WBC 5 part differential 3.0mL

ASSAY SHEET

LOT	MEK5D242L
	2024-05-05
OPEN	14DAYS

MEK-9100,9200
2D Barcode Entry



		RANGE	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200
WBC	10 ³ /μL	±0.6	3.1	3.1	3.1	3.07
NE%	%	±15.0	22.0	22.0	22.0	22.05
LY%	%	±15.0	56.9	56.9	56.9	61.63
MO%	%	±5.0	12.4	12.4	12.4	8.93
EO%	%	±5.0	4.9	4.9	4.9	5.69
BA%	%	±4.0	3.8	3.8	3.8	1.70
NE	10 ³ /μL	±0.5	0.6	0.6	0.6	0.69
LY	10 ³ /μL	±0.5	1.8	1.8	1.8	1.89
MO	10 ³ /μL	±0.2	0.4	0.4	0.4	0.27
EO	10 ³ /μL	±0.2	0.2	0.2	0.2	0.17
BA	10 ³ /μL	±0.2	0.1	0.1	0.1	0.05
RBC	10 ⁶ /μL	±0.15	2.59	2.57	2.60	2.54
HGB	g/dL	±0.4	6.6	6.7	6.5	6.70
HCT	%	±2.0	21.8	21.0	20.6	23.2
MCV	fL	±7.0	84.1	81.6	79.4	91.3
MCH	pg	±2.4	25.5	26.1	25.0	26.4
MCHC	g/dL	±3.0	30.3	31.9	31.6	28.9
RDW-CV *	%	±4.0	14.7	15.9	15.0	14.3
PLT	10 ³ /μL	±30	57	64	58	55.7
MPV	fL	±3.0	4.4	5.4	4.7	7.1

Manual Registration	
Lot	242
Level	Low
Control	MEK5D

* RDW-CV is displayed as RDW on MEK-7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901690N