


MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS242
	2024-03-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



		MEK -6400 -6410 -6420	MEK -6500 -6510	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305	MEK -6450	MEK -6550
WBC	10 ³ /μL	8.9	8.9	9.1	9.0	8.9	9.24	9.00	8.9	8.9
RBC	10 ⁶ /μL	4.70	4.70	4.70	4.67	4.70	4.75	4.69	4.70	4.70
HGB	g/dL	13.8	13.9	13.7	13.6	13.5	13.87	13.83	13.8	13.9
HCT	%	40.1	40.1	41.8	40.6	39.5	46.8	46.1	40.1	40.1
MCV	fL	85.3	85.3	88.9	86.9	84.1	98.6	98.4	85.3	85.3
MCH	pg	29.4	29.6	29.1	29.1	28.7	29.2	29.5	29.4	29.6
MCHC	g/dL	34.4	34.7	32.8	33.5	34.2	29.6	30.0	34.4	34.7
PLT	10 ³ /μL	258	258	248	251	250	265.1	271.0	258	258
MPV	fL	7.1	8.9	7.1	7.8	7.1	9.4	9.4	7.1	8.9
RDW-CV *	%	15.5	15.5	14.4	15.6	14.7	15.2	15.2	15.5	15.5

MEK	
	-7222
	-7300
	-8222
	-9100
	-9200
FS	160
FL	92
SD	102
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.



<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN



0694-901780L