


MK-3CN HEMATOLOGY CONTROL+CRP 2.5mL


ASSAY SHEET

LOT	MK3C245N
	2024-08-05
OPEN	14DAYS

MEK-1303 Barcode Entry

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		MEK-1303	RANGE
WBC	10 ³ /μL	7.64	± 0.80
RBC	10 ⁶ /μL	4.64	± 0.20
HGB	g/dL	13.86	± 0.50
HCT	%	45.0	± 3.0
MCV	fL	97.0	± 7.0
MCH	pg	29.9	± 2.4
MCHC	g/dL	30.8	± 3.0
PLT	10 ³ /μL	275.2	± 50.0
MPV	fL	8.3	± 2.0
RDW-CV	%	16.1	± 3.0
LY%	%	37.03	± 10.00
MO%	%	2.10	± 5.00
GR%	%	60.87	± 10.00
LY	10 ³ /μL	2.83	± 1.00
MO	10 ³ /μL	0.16	± 0.40
GR	10 ³ /μL	4.65	± 1.00
CRP	Precision	≤ 10.0%	
		(20.90-38.90 mg/L)	

Manual Registration	
Lot	245
Level	Normal
Control	MK-3C


<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION

 1-31-4 NISHIOCHIAI, SHINJUKU-KU
 TOKYO 161-8560 JAPAN

 0123
0694-901766F