


MK-RE3 HEMATOLOGY CONTROL
for reticulocyte **3.0mL**

ASSAY SHEET

LOT	MKRE245-3
	2024-07-05
OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	8.12	±3.88
RET	10 ⁶ /μL	0.2926	±0.1398

2D Barcode Entry



• If RBC value is under 3.00 [$\times 10^6/\mu\text{L}$], please refrain from using.
It may have deteriorated or hemolyzed, which affects the measurement accuracy.



<https://www.nihonkohden.com/spi/index.html>

Manual Registration					
Lot	245				
Level	3				
Expiration	7	M	5	D	2024 Y
Control	MKRE				



NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN