

MEK-CAL HEMATOLOGY CALIBRATOR 2.0mL

ASSAY SHEET

LOT	PLUS246
	2024-07-05
OPEN	7DAYS

MEK-9100,9200
2D Barcode Entry



MEK-1301,1302,1303,1305
Barcode Entry



		MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -7222	MEK -7300	MEK -8222	MEK -9100 -9200	MEK -1301 -1302 -1303 -1305
WBC	10 ³ /μL	9.1	9.1	9.2	9.1	9.0	9.28	9.08
RBC	10 ⁶ /μL	4.67	4.67	4.63	4.58	4.62	4.70	4.64
HGB	g/dL	13.8	13.9	13.6	13.6	13.5	13.80	13.82
HCT	%	40.3	40.3	41.6	40.3	39.2	46.7	46.1
MCV	fL	86.4	86.4	89.8	88.0	84.9	99.4	99.4
MCH	pg	29.6	29.8	29.4	29.7	29.2	29.4	29.8
MCHC	g/dL	34.2	34.5	32.7	33.7	34.4	29.6	30.0
PLT	10 ³ /μL	239	239	235	238	237	248.7	256.2
MPV	fL	6.0	8.4	6.0	7.3	6.0	8.7	8.7
RDW-CV *	%	15.5	15.5	14.6	15.9	14.9	15.5	15.4

MEK	-7222
	-7300
	-8222
	-9100
	-9200
FS	166
FL	89
SD	107
FS THR.	100

* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450,7222,8222.

<https://www.nihonkohden.com/spi/index.html>

NIHON KOHDEN CORPORATION
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN

MEK-CAL

PLUS246

2024-07-05

0123
0694-901780P