


MK-RE1 HEMATOLOGY CONTROL
for reticulocyte **3.0mL**

ASSAY SHEET

LOT	MKRE247-1
	2024-09-05
OPEN	14 DAYS

		MEK-9200	RANGE
RET%	%	1.56	±1.44
RET	10⁶/μL	0.0561	±0.0518

2D Barcode Entry



• If RBC value is under 3.00 [$\times 10^6/\mu\text{L}$], please refrain from using.
It may have deteriorated or hemolyzed, which affects the measurement accuracy.



<https://www.nihonkohden.com/spi/index.html>

Manual Registration						
Lot	247					
Level	1					
Expiration	9	M	5	D	2024	Y
Control	MKRE					



 **NIHON KOHDEN CORPORATION**
1-31-4 NISHIOCHIAI, SHINJUKU-KU
TOKYO 161-8560 JAPAN