


**MEK-3DH** HEMATOLOGY CONTROL 2.0mL

ASSAY SHEET

LOT	B248H
	2024-11-05
OPEN	14DAYS

**MEK-1301,1302,1305 Barcode Entry**

The barcodes may be read in any order. If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.

1



2



3



		RANGE	MEK -6400 -6410 -6420 -6450	MEK -6500 -6510 -6550	MEK -1301 -1302 -1305
WBC	10 <sup>3</sup> /μL	±2.2	21.1	21.1	20.61
RBC	10 <sup>6</sup> /μL	±0.25	5.91	5.91	5.94
HGB	g/dL	±0.7	19.2	19.1	19.37
HCT	%	±4.0	55.1	55.1	63.6
MCV	fL	±7.0	93.3	93.3	107.1
MCH	pg	±2.8	32.5	32.3	32.6
MCHC	g/dL	±3.0	34.8	34.7	30.5
PLT	10 <sup>3</sup> /μL	±70	530	530	548.9
MPV	fL	±3.0	5.7	8.1	8.7
RDW-CV *	%	±3.0	15.8	15.8	17.2
LY%	%	±10.0	16.7	16.7	13.18
MO%	%	±5.0	3.4	3.4	2.50
GR%	%	±10.0	79.9	79.9	84.32
LY	10 <sup>3</sup> /μL	±1.8	3.5	3.5	2.72
MO	10 <sup>3</sup> /μL	±0.7	0.7	0.7	0.52
GR	10 <sup>3</sup> /μL	±2.0	16.9	16.9	17.37

Manual Registration	
Lot	248
Level	High
Control	MEK-3D

\* RDW-CV is displayed as RDW on MEK-6400,6410,6420,6450.



<https://www.nihonkohden.com/spi/index.html>



**NIHON KOHDEN CORPORATION**  
1-31-4 NISHIOCHIAI, SHINJUKU-KU  
TOKYO 161-8560 JAPAN



0694-901686Q