


**MK-3CL** HEMATOLOGY CONTROL+CRP 2.5mL

**ASSAY SHEET**

<b>LOT</b>	<b>MK3C248L</b>
	<b>2024-11-05</b>
<b>OPEN</b>	<b>14DAYS</b>

**MEK-1303 Barcode Entry**

If the barcode is difficult to scan, place this assay sheet on a flat surface, keep the barcode reader and barcode vertical, and scan within a distance of 2cm.


1



2



3



		MEK-1303	RANGE
WBC	10 <sup>3</sup> /μL	2.01	± 0.60
RBC	10 <sup>6</sup> /μL	2.27	± 0.15
HGB	g/dL	6.02	± 0.40
HCT	%	20.2	± 2.0
MCV	fL	88.8	± 9.0
MCH	pg	26.5	± 2.4
MCHC	g/dL	29.8	± 3.0
PLT	10 <sup>3</sup> /μL	73.3	± 30.0
MPV	fL	8.6	± 3.0
RDW-CV	%	15.7	± 3.0
LY%	%	66.31	± 10.00
MO%	%	5.85	± 5.00
GR%	%	27.84	± 10.00
LY	10 <sup>3</sup> /μL	1.33	± 0.50
MO	10 <sup>3</sup> /μL	0.12	± 0.20
GR	10 <sup>3</sup> /μL	0.56	± 0.50
CRP	Precision	≤ 10.0%	
		(1.30-7.30 mg/L)	

Manual Registration	
Lot	<b>248</b>
Level	<b>Low</b>
Control	<b>MK-3C</b>


<https://www.nihonkohden.com/spi/index.html>

**NIHON KOHDEN CORPORATION**

 1-31-4 NISHIOCHIAI, SHINJUKU-KU  
 TOKYO 161-8560 JAPAN

**0694-901765G**